

## Measuring standard metabolic rates of insects by isothermal calorimetry

**Instruments to which this note applies:** I-Cal Flex

**Prepared by:** Lars Wadsö

**Target use:** ecology, insect metabolism, small animal activity, standard metabolic rate

### Introduction

The standard metabolic rate (SMR) of an animal is its lowest metabolic rate, its metabolic rate at rest. It is equivalent to basal metabolic rate (BMR), but while BMR usually is used for endothermic organisms that maintain a constant body temperature by their metabolic heat, SMR is a term used for poikilotherms, whose body temperature varies according to the ambient temperature. BMR has to be measured at thermal neutral conditions, while SMR can be measured at any temperature,

SMR can be reported in different ways, most commonly as thermal power (W), specific thermal power (W/g) or oxygen consumption rate. The SMRs based on thermal power (heat production rate) are most suitably measured by isothermal calorimetry.

For measurements of SMR to be of any value a few criteria have to be met:

- The measurement temperature needs to be stated.
- The measurement should be made while the animal is resting.
- The animal should be in good condition.
- The animal should not recently have been feeding.

SMRs are of interest for many reasons. One is that they are fundamental to theories of allometric relationships, i.e., how animal size and activity relate. The most famous of these relationships is Kleiber's law for endothermic animals, but there are also more recent and more general such models, for example Metabolic Theory of Ecology [1].

This application note concerns measurements of SMR on insects. As insects are primarily aerobic [2] and the enthalpy of respiration ( $\text{kJ/mol}_{\text{oxygen}}$ ) is independent of the type of substrate [3], calorimetric measurements of heat production rate of insects are equivalent to oxygen consumption measurements, but they do have a number of advantages compared to respirometry. The main advantage is that calorimeters continuously measure the rate as a function of time, so the resulting 'thermogram' will in many cases clearly show periods of rest and periods of activity, making it easy to determine the metabolic activity at rest.

### Materials and methods

Results from three different insects (Fig. 1) are presented in this application note:

1. Long-tailed silverfish (*Ctenolepisma longicaudata*) (length 18 mm, mass 60 mg) measured at 20 °C. Measurement was made with the insect in a plastic vial. Empty reference stack.
2. Larvae of a soldier beetle (family Cantharidae) (length 15 mm; mass 60 mg) measured at 20 °C.
3. Four larvae of the common cockchafer (*Melolontha melolontha*) (length 35 mm, mass 1.52-1.76 g) measured individually at 10 °C.



Fig. 1 – The insects investigated.

In all three cases baseline measurements were performed. This is important as the measured thermal powers are rather low. The three measurements were made in slightly different ways, with stages marked A, B and C, as seen in Fig. 2:

1. A: Baseline empty vial.  
B: Vial with insect.
2. A: Baseline empty vial.  
B: Soil and insect.  
C: Only soil.
3. A: Soil and insect.  
B: Only soil.

Soil was used in measurements 2 and 3 as both soldier beetle larvae and cockchafer larvae live in soil. It is important to keep the animals under conditions resembles what they are used to. This can in some cases be difficult, for example with flying insects or with social insects that are constantly active, trying to return to their colony [2].

### Results and discussion

Figure 2 shows the results, from which it is in all cases possible to calculate SMR, although in slightly different ways.

This animal is quite relaxed with




Fig. 2 – Result from measurement on soldier beetle larvae (A).

### 1. Long-tailed silverfish

Long tailed silverfish are found in warm places indoors in large parts of the world. The present measurement was made animal that had been kept at about 22 °C for several months. Its transfer to an empty vial does not seem to have give rise to stress as the animal was resting for the duration of the measurement (cf. Fig. 2B-C with activity periods) The SMR was found as the difference between the results from part B and part A (see Table 1).

### 2. Soldier beetle larvae

In this measurement a standard baseline was first run (A), then soil and the insect was added (B); finally the insect was removed and only the soil was measured (C). The soil was used to create good conditions for the larvae, but as soil has significant microbial activity, the thermal power from the soil had to be measured and used as a baseline above which the insect metabolism was found. From the results it is clear that periods of rest interleave periods of activity with higher thermal power. It is easy to construct a linear curve fit to the resting periods and compare this with a linear regression of the soil results, to find the SMR as the difference between these lines.

### 3. Cockchafer larvae

Here, measurements on four larvae were made simultaneously and it was found that one larvae had a quite different 'thermogram', with higher heat production, but not activity periods. The reason for this is not known. Also in this measurement is it easy to see the resting periods, and find the difference between the soil+insect (A) and soil (B) measurements.

Table 1 – The result

Sample	Thermal power	
	µW	µW/mg
1	53	0.88
2	145	2.4
3	58*	0.036*

\* The three lower results in Fig. 2C.

### Notes on the calorimetric measurement

#### Standard metabolic rate (SMR) measurements

It is not trivial to measure SMR [2] and one problem is that it is difficult to know if the animal was at rest during a respirometry measurement or not [4]. Researchers have tried to solve this issue, for example by using cameras to record activity periods, but it seems that isothermal calorimetry is a simpler way, as it is easy to pick out activity periods (Fig. 2B-C).

#### Conditions for animals during measurement

It is clear that many animals do not have their natural behavior when they are measured on in respirometers or calorimeters. However, it is possible to try to make the environment in, e.g., a calorimetric vial as good as possible, for example by using soil in the vial for soil living insects. Generally, the following points could be considered:

- The size of the animal in relation to the size of the vial.
- The activity and activity modes of the animal.
- The texture of the inside of the vial.
- If the round shape of normal vial is a problem (could induce circular motion).
- The oxygen and carbon dioxide levels during a measurement.

### References

1. Sibly, R.M., J.H. Brown, and A. Kodric-Brown, eds. *Metabolic ecology: a scaling approach*. 2012, Wiley-Blackwell.
2. Waters, J.S. and J.F. Harrison, *Insect metabolic rates*, in *Metabolic ecology - a scaling approach*, R.M. Sibly, J.H. Brown, and A. Kodric-Brown, Editors. 2012, Wiley-Blackwell. p. 198-211.
3. Wadsö, L. and L.D. Hansen, *Calorespirometry of terrestrial organisms and ecosystems*. *Methods*, 2015. 76(Online 28 October 2014): p. 11-19.
4. DeVries, Z.C. and A.G. Appel, *Standard metabolic rates of *Lepisma saccharina* and *Thermobia domestica*: Effects of temperature and mass*. *J. Insect Physiol.*, 2013. 59: p. 638-645.